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EN ORTHODONTIC PACIFIER THERMOMETER IN SILICONE

The orthodontic pacifier is soft and safe for the gums, easy to use, handy to carry and sold with a protection lid. Press on the key ON/OFF in order to switch on the Thermometer function, the last temperature reading measurement appears.

Temperature range: $Lo < 32,0^{\circ}C - 42,9^{\circ}C > Hi$

Accuracy: $\pm 0,1^{\circ}C$ for a temperature range between $35^{\circ}C$ to $42^{\circ}C$.

Temperature measurement timing : between 3 and 5 minutes

Device sealed, no dismountable. The battery cannot be replaced.

Beep at the beginning and end of the temperature measurement. Sound alarm in case of fever (above $37,8^{\circ}C$).

The indicator "Err" appears in case of a dysfunction of the system.

Reference standards Clinical Thermometers EN12470-3 : 2000

Electromagnetic compatibility EMC EN60601-1-2 and biocompatibility test in accordance with Annex 1 of the MDD

1 year warranty under normal use

CLEANING AND DISINFECTION

Please make sure to sterilize the pacifier after each use. You can clean it either with hot and soapy water or wipe it with a cotton cloth soaked in 70% alcohol. DO NOT BOIL STEAM OR WASH IN THE DISHWASHER. It would damage your thermometer and cancel the warranty.

DO NOT LEAVE WITHIN THE REACH OF AN UNATTENDED CHILD

HOW TO USE

Maintain the key ON/OFF pushed until seeing the indication " $188,8^{\circ}E/M$ " on the LCD screen. Release it, the screen displays $Lo^{\circ}C/M$, then $37^{\circ}C$ and finally $Lo^{\circ}C$.

At the beep, place the pacifier into the child's mouth and let him/her suck the thermometer pacifier for at least 3 minutes. A new beep indicates the end of the temperature measurement. Withdraw the pacifier and read the temperature measurement on the screen. Automatic shutdown ($\pm 9min$).

Do not discard your used batteries with your domestic waste. Bring them to a collection organization. In case of an extended non use, please remove the batteries from the device in order to avoid any liquid leak that could damage your device.