



Pos. Art.Nr. Anz. Bezeichnung

- | | | | |
|----|----------|----|---|
| 1 | DR200-KC | 1 | Potassium chloride, 250 g |
| 2 | DR201-1C | 1 | Columbit |
| 3 | DR250-1A | 1 | Radiation absorption plates, set - consisting of: |
| | | 5 | Aluminium |
| | | 5 | Iron |
| | | 5 | Acryl glass |
| | | 10 | Lead |
| 4 | DR212-1H | 1 | Absorption plate mount, magnetic |
| 5 | DR201-1R | 1 | Radioactive preparation mount, magnetic |
| 6 | C3551-2T | 1 | Test tube, graduated |
| 7 | DM115-1A | 1 | Lead (tare) shot, 250 g |
| 8 | DR213-1A | 1 | Adapter for deflection in radioactive substances |
| 9 | DE407-1A | 1 | Magnets, button type, pair of |
| 10 | C6008-1B | 1 | Sample box with lid, plastics, 80 ml, 50 x 50 x 40 mm |
| 11 | C7418-2A | 1 | Knife for laboratory use, steel |
| 12 | C7415-2Z | 1 | Crucible tongs, stainless steel, L = 200 mm |

Aufbewahrung:

- | | | | |
|----|----------|---|------------------------------------|
| 13 | P7906-1R | 1 | Box-insert Radioactivity, plastics |
| 14 | P7806-1K | 1 | Storage box II small, with cover |

Plan for storage of all components